

产品名称: TSNAX Mouse Monoclonal Antibody
产品货号: AMM83622



产品概述 (Summary)

产品名称 (Production Name)	TSNAX Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TSNAX
别名 (Alternative Names)	TSNAX
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q99598.Purified recombinant fragment of human TSNAX

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:250
蛋白分子量 (Molecular Weight)	Calculated MW: 33 kDa; Observed MW: 33 kDa


研究背景 (Background)

产品名称: TSNAX Mouse Monoclonal Antibody
产品货号: AMM83622




研究领域 (Research Area)

图片 (Image Data)



Lane 1: 15 μ g of MCF7 lysate; Lane 2: 15 μ g of HeLa lysate; Lane 3: 15 μ g of HEK-293 lysate; Lane 4: 15 μ g of SK-OV-3 lysate;

Result: AMM83622 can detect TSNAX by Western blotting.



IHC-P analysis of testis tissue by TSNAX antibody (AMM83622)

注意事项 (Note)

For research use only .