

产品名称: ASPM Mouse Monoclonal Antibody
产品货号: AMM83623



产品概述 (Summary)

产品名称 (Production Name)	ASPM Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ASPM
别名 (Alternative Names)	ASPM
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q8IZT6.Purified recombinant fragment of human ASPM

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:250
蛋白分子量 (Molecular Weight)	382kDa


研究背景 (Background)

产品名称: ASPM Mouse Monoclonal Antibody
产品货号: AMM83623




研究领域 (Research Area)

图片 (Image Data)



Western blot analysis using ASPM mouse mAb against HEK-293 lysate



IHC-P analysis of adrenal gland tissue by ASPM antibody (AMM83623)

注意事项 (Note)

For research use only .