

产品名称: IL18 Mouse Monoclonal Antibody
产品货号: AMM83625



产品概述 (Summary)

产品名称 (Production Name)	IL18 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IL18
别名 (Alternative Names)	IL18
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q14116.Purified recombinant fragment of human IL18

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:100-1:200,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 22 kDa; Observed MW: 17 kDa


研究背景 (Background)

产品名称: IL18 Mouse Monoclonal Antibody
产品货号: AMM83625




研究领域 (Research Area)


图片 (Image Data)



Western blot analysis using IL18 mouse mAb against U-2 OS lysate



Lane 1: U-2 OS lysate; Lane 2: IL18 immunoprecipitated from U-2 OS lysate by AMM83625; Lane 3: The same as Lane 2 but use IgG isotype control antibody; Result: AMM83625 can immunoprecipitate IL18;



IHC-P analysis of tonsil tissue by IL18 antibody (AMM83625)

注意事项 (Note)

产品名称: IL18 Mouse Monoclonal Antibody
产品货号: AMM83625



For research use only .