

产品名称: HIP1 Mouse Monoclonal Antibody
产品货号: AMM83626



产品概述 (Summary)

产品名称 (Production Name)	HIP1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HIP1
别名 (Alternative Names)	HIP1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O00291.Purified recombinant fragment of human HIP1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:100-1:200,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	116.2kDa


研究背景 (Background)

产品名称: HIP1 Mouse Monoclonal Antibody
产品货号: AMM83626




研究领域 (Research Area)

图片 (Image Data)




Lane 1: 15 µg of A549 lysate; Lane 2: 15 µg of HCT116 lysate; Lane 3: 15 µg of HEK-293 lysate; Lane 4: 15 µg of HeLa lysate;

Result: AMM83626 can detect HIP1 by Western blotting.



Lane 1: A549 lysate; Lane 2: HIP1 immunoprecipitated from A549 lysate by AMM83626; Lane 3: The same as Lane 2 but use IgG isotype control antibody; Result: AMM83626 can immunoprecipitate HIP1;



IHC-P analysis of kidney tissue by HIP1 antibody (AMM83626)

产品名称: HIP1 Mouse Monoclonal Antibody
产品货号: AMM83626



注意事项 (Note)

For research use only .