

产品名称: GMPR2 Mouse Monoclonal Antibody
产品货号: AMM83629



产品概述 (Summary)

产品名称 (Production Name)	GMPR2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GMPR2
别名 (Alternative Names)	GMPR2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9P2T1.Purified recombinant fragment of human GMPR2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	37.9kDa


研究背景 (Background)

产品名称: GMPR2 Mouse Monoclonal Antibody
产品货号: AMM83629




研究领域 (Research Area)

图片 (Image Data)



Western blot analysis using GMPR2 mouse mAb against Hep G2 lysate



Lane 1: GMPR2 immunoprecipitated from Hep G2 lysate by AMM83629; Lane 2: The same as Lane 2 but use IgG isotype control antibody; Result: AMM83629 can immunoprecipitate GMPR2.

注意事项 (Note)

For research use only .