

产品名称: ECI2 Mouse Monoclonal Antibody
产品货号: AMM83663



产品概述 (Summary)

产品名称 (Production Name)	ECI2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ECI2
别名 (Alternative Names)	ECI2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O75521.Purified recombinant fragment of human ECI2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:100-1:200,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	40kDa


研究背景 (Background)

产品名称: ECI2 Mouse Monoclonal Antibody
产品货号: AMM83663




研究领域 (Research Area)


图片 (Image Data)



Western blot analysis using ECI2 mouse mAb against Raji lysate



Lane 1: Raji lysate ; Lane 2: ECI2 immunoprecipitated from Raji lysate by AMM83663 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83663 can immunoprecipitate ECI2;



IHC-P analysis of liver tissue by ECI2 antibody (AMM83663)

注意事项 (Note)

产品名称: ECI2 Mouse Monoclonal Antibody
产品货号: AMM83663



For research use only .