

产品名称: LAT2 Mouse Monoclonal Antibody
产品货号: AMM83664



产品概述 (Summary)

产品名称 (Production Name)	LAT2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LAT2
别名 (Alternative Names)	LAT2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9GZY6.Purified recombinant fragment of human LAT2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:20-1:100,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	26.6kDa


研究背景 (Background)

产品名称: LAT2 Mouse Monoclonal Antibody
产品货号: AMM83664




研究领域 (Research Area)


图片 (Image Data)



Western blot analysis using LAT2 mouse mAb against U-937 lysate



Lane 1: U-937 lysate ; Lane 2: LAT2 immunoprecipitated from U-937 lysate by AMM83664 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83664 can immunoprecipitate LAT2;



IHC-P analysis of tonsil tissue by LAT2 antibody (AMM83664)

注意事项 (Note)

产品名称: LAT2 Mouse Monoclonal Antibody
产品货号: AMM83664



For research use only .