

产品名称: IRAG2 Mouse Monoclonal Antibody
产品货号: AMM83665



产品概述 (Summary)

产品名称 (Production Name)	IRAG2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IRAG2
别名 (Alternative Names)	IRAG2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q12912.Purified recombinant fragment of human IRAG2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:20-1:100,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:;Observed MW:61kD


研究背景 (Background)

产品名称: IRAG2 Mouse Monoclonal Antibody
产品货号: AMM83665




研究领域 (Research Area)


图片 (Image Data)



Western blot analysis using IRAG2 mouse mAb against Ramos lysate



Lane 1: Ramos lysate ; Lane 2: IRAG2 immunoprecipitated from Ramos lysate by AMM83665 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83665 can immunoprecipitate IRAG2;



IHC-P analysis of tonsil tissue by IRAG2 antibody (AMM83665)

注意事项 (Note)

产品名称: IRAG2 Mouse Monoclonal Antibody
产品货号: AMM83665



For research use only .