

产品名称: Acetyl-Alpha Tubulin (Lys40) Mouse Monoclonal Antibody
产品货号: AMM84825



产品概述 (Summary)

产品名称 (Production Name)	Acetyl-Alpha Tubulin (Lys40) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Acetylated
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.5%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Acetyl-Alpha Tubulin (Lys40)
别名 (Alternative Names)	alpha tubulin 2; B alpha 1; H2 alpha; Hum a tub 1; TUBA1; TUBA1A; Tubulin alpha 1chain
基因 ID (Gene ID)	10376.0
蛋白 ID (SwissProt ID)	P68363.Synthetic peptide conjugated to KLH.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 50 kDa; Observed MW: 50 kDa

产品名称: Acetyl-Alpha Tubulin (Lys40) Mouse Monoclonal Antibody
产品货号: AMM84825




研究背景 (Background)


TUBA1A Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable site on the beta chain and one at a non-exchangeable site on the alpha-chain.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of Acetyl- α -Tubulin (Lys40) in rat spleen, Raw246.7 lysates using Acetyl- α -Tubulin (Lys40) antibody.



Western blot analysis of AcetylAlpha Tubulin (Lys40) in 293, C6 lysates using Acetyl- α -Tubulin (Lys40) antibody.



产品名称: Acetyl-Alpha Tubulin (Lys40) Mouse Monoclonal Antibody
产品货号: AMM84825



Immunohistochemistry analysis of paraffin-embedded Human Brain Tissue using α -tubulin(Acetyl Lys4) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .