

产品名称: CD23 Mouse Monoclonal Antibody
产品货号: AMM84934



产品概述 (Summary)

产品名称 (Production Name)	CD23 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC, ICC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD23 FCER2; CD23A; CLEC4J; FCE2; IGEBF; Low affinity immunoglobulin epsilon Fc receptor; BLAST-2; C-type lectin domain family 4 member J; Fc-epsilon-RII; Immunoglobulin E-binding factor; Lymphocyte IgE receptor; CD23
别名 (Alternative Names)	
基因 ID (Gene ID)	2208.0
蛋白 ID (SwissProt ID)	P06734.Synthetic Peptide of CD23

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:100, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	/

产品名称: CD23 Mouse Monoclonal Antibody
产品货号: AMM84934




研究背景 (Background)


This receptor has essential roles in the regulation of IgE production and in the differentiation of B-cells (it is a B-cell-specific antigen).

研究领域 (Research Area)


图片 (Image Data)



Immunofluorescence analysis of CD23 in Human stomach using CD23 antibody(red),and DAPI (blue).



Immunofluorescence analysis of CD23 in rat lung tissue using CD23 antibody(1E9)(red),and DAPI (blue).




Immunohistochemistry analysis of paraffin-embedded Human lung tissue using CD23 antibody. High-pressure and


产品名称: CD23 Mouse Monoclonal Antibody
产品货号: AMM84934




temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.



Immunohistochemistry analysis of paraffin-embedded rat kidney tissue using CD23 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.



Immunohistochemistry analysis of paraffin-embedded mouse lung tissue using CD23 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.



Immunohistochemistry analysis of paraffin-embedded Human tonsil tissue using CD23 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .