

产品名称: MYH1 Mouse Monoclonal Antibody
产品货号: AMM84935



产品概述 (Summary)

产品名称 (Production Name)	MYH1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ICC
种属反应性 (Reactivity)	Human,Mouse,Rat,Fruitfly,Nematode

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.5%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MYH1 fast skeletal muscle; MYH3; HEMHC; Myosin-1; Myosin heavy chain 2x; MyHC-2x; Myosin heavy chain IIx/d; MyHC-IIx/d; Myosin heavy chain; skeletal muscle; adult 1
别名 (Alternative Names)	
基因 ID (Gene ID)	4619.0
蛋白 ID (SwissProt ID)	P12882.Synthetic Peptide of Myosin Heavy Chain

产品应用 (Application)

稀释比 (Dilution Ratio)	ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	/

产品名称: MYH1 Mouse Monoclonal Antibody
产品货号: AMM84935




研究背景 (Background)

Muscle contraction.

研究领域 (Research Area)

图片 (Image Data)



Immunofluorescence analysis of MYH1 in Immunofluorescence Staining of nematode tissue using MYH1 antibody.

注意事项 (Note)

For research use only .