

产品名称: IL-2 Receptor alpha Mouse Monoclonal Antibody
产品货号: AMM84939



产品概述 (Summary)

产品名称 (Production Name)	IL-2 Receptor alpha Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IL-2 Receptor alpha
别名 (Alternative Names)	IL2RA; Interleukin-2 receptor subunit alpha; IL-2 receptor subunit alpha; IL-2-RA; IL-2R subunit alpha; IL2-RA; TAC antigen; p55; CD25
基因 ID (Gene ID)	3559.0
蛋白 ID (SwissProt ID)	P01589.Synthetic Peptide of CD25

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	/

研究背景 (Background)

产品名称: IL-2 Receptor alpha Mouse Monoclonal Antibody
产品货号: AMM84939




Receptor for interleukin-2.


研究领域 (Research Area)

TGF-beta signaling pathway, PI3K-Akt signaling pathway, Jak-STAT signaling pathway

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded Human liver tissue using IL-2 Receptor alpha antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.



Immunohistochemical analysis of paraffin-embedded Human tonsils using IL-2 Receptor alpha antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .