

产品名称: Progesterone Receptor Mouse Monoclonal Antibody
产品货号: AMM84940



产品概述 (Summary)

产品名称 (Production Name)	Progesterone Receptor Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC, ICC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Progesterone Receptor
别名 (Alternative Names)	PGR; NR3C3; Progesterone receptor; PR; Nuclear receptor subfamily 3 group C member 3
基因 ID (Gene ID)	5241.0
蛋白 ID (SwissProt ID)	P06401.Synthetic Peptide of PR

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:100, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	/

产品名称: Progesterone Receptor Mouse Monoclonal Antibody
产品货号: AMM84940




研究背景 (Background)


Progesterone plays a central role in the reproductive events associated with the establishment and maintenance of pregnancy. Progesterone receptor, a member of the steroid receptor superfamily, mediates the physiologic effects of progesterone.

研究领域 (Research Area)

图片 (Image Data)




Immunofluorescence analysis of Progesterone Receptor in Human appendix tissue using Progesterone Receptor antibody(Z15)(red),and DAPI (blue).




Immunofluorescence analysis of Progesterone Receptor in mouse liver using PR antibody(Z15)(red),and DAPI (blue).


产品名称: Progesterone Receptor Mouse Monoclonal Antibody
产品货号: AMM84940



Immunofluorescence analysis of Progesterone Receptor in rat heart using Progesterone Receptor antibody(red), and DAPI (blue).




Immunohistochemistry analysis of paraffin-embedded Human uterus tissue using Progesterone Receptor antibody.High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.Negative control was used by secondary antibody only.



Immunohistochemical analysis of paraffin-embedded Human tonsils using Progesterone Receptor antibody.High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.Negative control was used by secondary antibody only.

产品名称: Progesterone Receptor Mouse Monoclonal Antibody
产品货号: AMM84940



Immunohistochemistry analysis of paraffin-embedded mouse testis tissue using PR antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.

注意事项 (Note)

For research use only .