

产品名称: alpha Actinin Mouse Monoclonal Antibody  
产品货号: AMM84946



## 产品概述 (Summary)

|                        |                                         |
|------------------------|-----------------------------------------|
| 产品名称 (Production Name) | alpha Actinin Mouse Monoclonal Antibody |
| 描述 (Description)       | Mouse monoclonal Antibody               |
| 宿主 (Host)              | Mouse                                   |
| 应用 (Application)       | IHC                                     |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat                         |

## 产品性能 (Performance)

|                     |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                |
| 修饰 (Modification)   | Unmodified                                                                                  |
| 同种型 (Isotype)       | Mouse IgG1                                                                                  |
| 克隆 (Clonality)      | Monoclonal                                                                                  |
| 形式 (Form)           | Liquid                                                                                      |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.    |
| 储存溶液 (Buffer)       | Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol. |
| 纯化方式 (Purification) | Affinity Purification                                                                       |

## 免疫原信息 (Immunogen)

|                        |                                                            |
|------------------------|------------------------------------------------------------|
| 基因名 (Gene Name)        | alpha Actinin                                              |
| 别名 (Alternative Names) | actinin 1 smooth muscle; actinin alpha 1; ACTN 1; FLJ40884 |
| 基因 ID (Gene ID)        | 87.0                                                       |
| 蛋白 ID (SwissProt ID)   | P12814. Purified recombinant protein expressed in E.coli.  |

## 产品应用 (Application)

|                          |                |
|--------------------------|----------------|
| 稀释比 (Dilution Ratio)     | IHC 1:50-1:100 |
| 蛋白分子量 (Molecular Weight) | /              |

## 研究背景 (Background)

F-actin cross-linking protein which is thought to anchor actin to a variety of intracellular structures. This is a bundling


**产品名称: alpha Actinin Mouse Monoclonal Antibody**  
**产品货号: AMM84946**




protein.

## **研究领域 (Research Area)**

## **图片 (Image Data)**



Immunohistochemistry analysis of paraffin-embedded Human Breast Carcinoma Tissue using alpha Actinin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using alpha Actinin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

## **注意事项 (Note)**

For research use only .