

产品名称: Cytokeratin 6 Mouse Monoclonal Antibody
产品货号: AMM84964



产品概述 (Summary)

产品名称 (Production Name)	Cytokeratin 6 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Cytokeratin 6 CK 6A; CK 6B; CK 6C; CK 6D; CK 6E; CK-6C; CK-6E; K6a keratin; K6b keratin; K6C; K6c keratin; K6d keratin; K6e keratin; Keratin; KRT6A; KRT6B; KRT6C; KRT6D; KRT6E
别名 (Alternative Names)	
基因 ID (Gene ID)	3853/3854/286887
蛋白 ID (SwissProt ID)	P02538/P04259/P48668.Synthetic Peptide of Cytokeratin 6

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IHC 1:50-1:100, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 60 kDa; Observed MW: 56-60 kDa

产品名称: Cytokeratin 6 Mouse Monoclonal Antibody
产品货号: AMM84964




研究背景 (Background)


There are at least six isoforms of human type II keratin-6 (K6). There are two types of cytoskeletal and microfibrillar keratin, I (acidic) and II (neutral to basic) (40-55 and 56-70 kDa, respectively).

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of Cytokeratin 6 in Hela, C6, 3T3, K562, HepG2 lysates using Cytokeratin 6 antibody



Immunofluorescence analysis of Cytokeratin 6 in Human breastcancer tissue using Cytokeratin 6 antibody(red),and DAPI (blue).



产品名称: Cytokeratin 6 Mouse Monoclonal Antibody
产品货号: AMM84964



Immunohistochemistry analysis of paraffin-embedded Human breastcancer tissue using Cytokeratin 6 antibody.High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.Negative control was used by secondary antibody only.

注意事项 (Note)

For research use only .