

产品名称: Rabbit IgG Light Chain Mouse Monoclonal Antibody
产品货号: AMM84973



产品概述 (Summary)

产品名称 (Production Name)	Rabbit IgG Light Chain Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Rabbit

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Rabbit IgG Light Chain
别名 (Alternative Names)	Antibody; Ig gamma chain C region; Ig kappa-b4 chain C region; Ig kappa-b5 chain C region; Ig kappa-b9 chain C region; Ig lambda chain C region; K-BAS.
基因 ID (Gene ID)	-
蛋白 ID (SwissProt ID)	-Recombinant Protein of Rabbit IgG (Light Chain Specific)

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 25 kDa; Observed MW: 25 kDa

产品名称: Rabbit IgG Light Chain Mouse Monoclonal Antibody
产品货号: AMM84973




研究背景 (Background)

Constant region of immunoglobulin light chains. Immunoglobulins, also known as antibodies, are membrane-bound or secreted glycoproteins produced by B lymphocytes.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Reduced Rabbit IgG, H+L using Rabbit IgG Light Chain antibody.

注意事项 (Note)

For research use only .