

产品名称: Prealbumin Mouse Monoclonal Antibody
产品货号: AMM84983



产品概述 (Summary)

产品名称 (Production Name)	Prealbumin Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Prealbumin
别名 (Alternative Names)	Prealbumin; TBPA; Amyloidosis I; Transthyretin
基因 ID (Gene ID)	7276.0
蛋白 ID (SwissProt ID)	P02766. Purified recombinant protein expressed in E.coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 16 kDa; Observed MW: 16 kDa

研究背景 (Background)


Thyroid hormone-binding protein. Probably transports thyroxine from the bloodstream to the brain.

产品名称: Prealbumin Mouse Monoclonal Antibody
产品货号: AMM84983




研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Prealbumin in Hela, Jurkat, mouse brain, MOLT4, Human serum lysates using Prealbumin antibody



Immunohistochemistry analysis of paraffin-embedded Human Liver Carcinoma Tissue using Prealbumin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using Prealbumin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: Prealbumin Mouse Monoclonal Antibody
产品货号: AMM84983



注意事项 (Note)

For research use only .