

产品名称: eIF2A Mouse Monoclonal Antibody
产品货号: AMM84996



产品概述 (Summary)

产品名称 (Production Name)	eIF2A Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB, ICC, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	eIF2A
别名 (Alternative Names)	EIF2S1; EIF2A; Eukaryotic translation initiation factor 2 subunit 1; Eukaryotic translation initiation factor 2 subunit alpha; eIF-2-alpha; eIF-2A; eIF-2alpha
基因 ID (Gene ID)	1965.0
蛋白 ID (SwissProt ID)	P05198. Purified recombinant human eIF2α protein fragments expressed in E.coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:50-1:200, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 36 kDa; Observed MW: 36 kDa

产品名称: eIF2A Mouse Monoclonal Antibody
产品货号: AMM84996




研究背景 (Background)


eIF2A a translation initiation factor that functions in the early steps of protein synthesis by forming a ternary complex with GTP and initiator tRNA. This complex binds to a 40s ribosomal subunit, followed by mRNA binding to form a 43S preinitiation complex.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of eIF2α in K562, COS7, 3T3, C2C12, MCF-7, C6 and HeLa lysates using eIF2α antibody.



Immunoprecipitation analysis of eIF2A in HeLa lysates using eIF2α antibody.



产品名称: eIF2A Mouse Monoclonal Antibody
产品货号: AMM84996



Immunofluorescence analysis of eIF2A in Hela using eIF2 α antibody.

注意事项 (Note)

For research use only .