

产品名称: Wnt5b Mouse Monoclonal Antibody
产品货号: AMM85002



产品概述 (Summary)

产品名称 (Production Name)	Wnt5b Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Transfected

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Wnt5b
别名 (Alternative Names)	Wnt5b
基因 ID (Gene ID)	81029.0
蛋白 ID (SwissProt ID)	Q9H1J7. Purified recombinant human Wnt5b protein fragments expressed in E.coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 40 kDa; Observed MW: 45 kDa

研究背景 (Background)

产品名称: Wnt5b Mouse Monoclonal Antibody
产品货号: AMM85002




Ligand for members of the frizzled family of seven transmembrane receptors. Probable developmental protein. May be a signaling molecule which affects the development of discrete regions of tissues. Is likely to signal over only few cell diameters.

研究领域 (Research Area)

Wnt signaling pathway

图片 (Image Data)



Western blot analysis of Wnt5b in CHO-K1 lysates (B) and CHO-K1 transfected by Wnt5b (A) lysates using Wnt5b antibody.

注意事项 (Note)

For research use only .