

产品名称: Wnt3 Mouse Monoclonal Antibody
产品货号: AMM85003



产品概述 (Summary)

产品名称 (Production Name)	Wnt3 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Transfected

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Wnt3
别名 (Alternative Names)	Proto-oncogene Wnt-3; INT4; TETAMS
基因 ID (Gene ID)	7473.0
蛋白 ID (SwissProt ID)	P56703. Purified recombinant human Wnt3 protein fragments expressed in E.coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 40 kDa; Observed MW: 40 kDa

研究背景 (Background)

产品名称: Wnt3 Mouse Monoclonal Antibody
产品货号: AMM85003




Ligand for members of the frizzled family of seven transmembrane receptors. Wnt-3 and Wnt-3a play distinct roles in cell-cell signaling during morphogenesis of the developing neural tube.

研究领域 (Research Area)

Wnt signaling pathway

图片 (Image Data)



Western blot analysis of Wnt3 in CHO-K1 lysates (B) and CHO-K1 transfected by Wnt3 (A) lysates using Wnt3 antibody.

注意事项 (Note)

For research use only .