

产品名称: RPS6 Mouse Monoclonal Antibody
产品货号: AMM85008



产品概述 (Summary)

产品名称 (Production Name)	RPS6 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB, ICC, IP
种属反应性 (Reactivity)	Human, Mouse, Rat, Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RPS6
别名 (Alternative Names)	RPS6; OK/SW-cl.2; 40S ribosomal protein S6; Phosphoprotein NP33
基因 ID (Gene ID)	6194.0
蛋白 ID (SwissProt ID)	P62753. Recombinant full length protein corresponding to human S6 Ribosomal Protein.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:50-1:200, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 29 kDa; Observed MW: 32 kDa

研究背景 (Background)

产品名称: RPS6 Mouse Monoclonal Antibody
产品货号: AMM85008




May play an important role in controlling cell growth and proliferation through the selective translation of particular classes of mRNA.


研究领域 (Research Area)

PI3K-Akt signaling pathway


图片 (Image Data)



Western blot analysis of S6 Ribosomal Protein in Ramos, Jurkat, C6, 3T3 and Hela lysates using S6 Ribosomal Protein antibody.




Western blot analysis of RPS6 in CHO-K1, Raji, Raw264.7 and COS7 lysates using RPS6 antibody



Immunoprecipitation analysis of RPS6 in Hela lysates using S6 Ribosomal Protein antibody.

产品名称: RPS6 Mouse Monoclonal Antibody
产品货号: AMM85008



Immunocytochemistry analysis of RPS6 in HeLa using S6 Ribosomal Protein antibody.

注意事项 (Note)

For research use only .