

产品名称: SSB Mouse Monoclonal Antibody
产品货号: AMM85010



产品概述 (Summary)

产品名称 (Production Name)	SSB Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SSB
别名 (Alternative Names)	La; LARP3; La/SSB
基因 ID (Gene ID)	6741.0
蛋白 ID (SwissProt ID)	P05455. Purified recombinant human La Antigen protein fragments expressed in E.coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 47 kDa; Observed MW: 50 kDa

研究背景 (Background)


产品名称: SSB Mouse Monoclonal Antibody
产品货号: AMM85010



Binds to the 3' poly(U) termini of nascent RNA polymerase III transcripts, protecting them from exonuclease digestion and facilitating their folding and maturation.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of La Antigen in MCF-7, A549, Jurkat and HeLa lysates using La Antigen antibody.

注意事项 (Note)

For research use only .