

产品名称: SOD2 Mouse Monoclonal Antibody
产品货号: AMM85014



产品概述 (Summary)

产品名称 (Production Name)	SOD2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SOD2
别名 (Alternative Names)	IPOB; IPO-B; MNSOD; MVCD6; Mn-SOD; SOD2
基因 ID (Gene ID)	6648.0
蛋白 ID (SwissProt ID)	P04179. Synthetic peptide corresponding to human SOD2 protein

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 25 kDa; Observed MW: 22 kDa

研究背景 (Background)

SOD-2 is a homotetrameric manganese enzyme (also known as MnSOD) that functions in the mitochondrion. ROS are


产品名称: SOD2 Mouse Monoclonal Antibody
产品货号: AMM85014



implicated in a wide range of degenerative processes, including Alzheimer' s disease, Parkinson' s disease and ischemic heart disease. Homozygous mutant mice, which lack SOD-2, exhibit dilated cardiomyopathy, accumulation of lipid in liver and skeletal muscle, metabolic acidosis, oxidative DNA damage and respiratory chain deficiencies in heart and skeletal muscle.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of SOD2 in HeLa, THP-1, A549 Ramos and K562 lysates using SOD2 antibody.

注意事项 (Note)

For research use only .