

产品名称: Cytokeratin 19 Mouse Monoclonal Antibody
产品货号: AMM85015



产品概述 (Summary)

产品名称 (Production Name)	Cytokeratin 19 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Cytokeratin 19
别名 (Alternative Names)	KRT19; Keratin; type I cytoskeletal 19; Cytokeratin-19; CK-19; Keratin-19; K19
基因 ID (Gene ID)	3880.0
蛋白 ID (SwissProt ID)	P08727.Synthetic Peptide of CK19

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 44 kDa; Observed MW: 44 kDa

研究背景 (Background)

Involved in the organization of myofibers. Together with KRT8, helps to link the contractile apparatus to dystrophin at the


产品名称: Cytokeratin 19 Mouse Monoclonal Antibody
产品货号: AMM85015




costameres of striated muscle.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Cytokeratin 19 in PC-3, SKOV-3 and MCF-7 lysates using Cytokeratin 19 antibody



Immunohistochemistry analysis of paraffin-embedded staining of Human colon cancer tissue using Cytokeratin 19 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .