

产品名称: HP1 gamma Mouse Monoclonal Antibody
产品货号: AMM85033



产品概述 (Summary)

产品名称 (Production Name)	HP1 gamma Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.5%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HP1 gamma
别名 (Alternative Names)	CBX3; Chromobox protein homolog 3; HECH; Heterochromatin protein 1 homolog gamma; HP1 gamma; Modifier 2 protein
基因 ID (Gene ID)	11335.0
蛋白 ID (SwissProt ID)	Q13185.Purified recombinant protein expressed in E.coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 21 kDa; Observed MW: 24 kDa

研究背景 (Background)


产品名称: HP1 gamma Mouse Monoclonal Antibody
产品货号: AMM85033




Recognizes and binds histone H3 tails methylated at 'Lys-9', leading to epigenetic repression. May contribute to the association of the heterochromatin with the inner nuclear membrane through its interaction with lamin B receptor (LBR). Involved in the formation of functional kinetochore through interaction with MIS12 complex proteins.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of HP1 gamma in HeLa lysates , 3T3 lysates , PC-12 lysates using HP1 gamma antibody



Western blot analysis of HP1 gamma in MCF-7, HepG2, Raw264.7, HeLa lysates using HP1 gamma antibody.




Immunohistochemistry analysis of paraffin-embedded Human Colon Carcinoma Tissue using HP1 gamma antibody.High-

产品名称: HP1 gamma Mouse Monoclonal Antibody
产品货号: AMM85033



pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using HP1 gamma antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .