

产品名称: Hsp40 Mouse Monoclonal Antibody
产品货号: AMM85034



产品概述 (Summary)

产品名称 (Production Name)	Hsp40 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Hsp40 DNAJB1; DNAJ1; HDJ1; HSPF1; DnaJ homolog subfamily B member 1; DnaJ protein homolog 1; Heat shock 40 kDa protein 1; HSP40; Heat shock protein 40; Human DnaJ protein 1; hDj-1
别名 (Alternative Names)	
基因 ID (Gene ID)	3337.0
蛋白 ID (SwissProt ID)	P25685.Recombinant Protein of Hsp40

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 38 kDa; Observed MW: 38 kDa

产品名称: Hsp40 Mouse Monoclonal Antibody
产品货号: AMM85034




研究背景 (Background)

Interacts with HSP70 and can stimulate its ATPase activity. Stimulates the association between HSC70 and HIP. Negatively regulates heat shock-induced HSF1 transcriptional activity during the attenuation and recovery phase period of the heat shock response (PubMed:9499401). Stimulates ATP hydrolysis and the folding of unfolded proteins mediated by HSPA1A/B (in vitro).

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of HSP40 in Jurkat lysates using HSP40 antibody.

注意事项 (Note)

For research use only .