

产品名称: AMPK beta 1 Mouse Monoclonal Antibody
产品货号: AMM85035



产品概述 (Summary)

产品名称 (Production Name)	AMPK beta 1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AMPK beta 1
别名 (Alternative Names)	5' '-AMP-activated protein kinase subunit beta-1; AMP-activated; noncatalytic; beta-1; AMPK; AMPK beta 1 chain; AMPK subunit beta-1; AMPK-BETA-1; AMPKb; HAMPKb; PRKAB1
基因 ID (Gene ID)	5564.0
蛋白 ID (SwissProt ID)	Q9Y478.Recombinant Protein of AMPK β 1.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 30 kDa; Observed MW: 38 kDa

产品名称: AMPK beta 1 Mouse Monoclonal Antibody
产品货号: AMM85035




研究背景 (Background)

AMP-activated protein kinase (AMPK) is highly conserved from yeast to plants and animals and plays a key role in the regulation of energy homeostasis. AMPK is a heterotrimeric complex composed of a catalytic α subunit and regulatory β and γ subunits, each of which is encoded by two or three distinct genes ($\alpha 1, 2$; $\beta 1, 2$; $\gamma 1, 2, 3$).

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of AMPK β 1 in HeLa, COS7 lysates using AMPK beta 1 antibody.

注意事项 (Note)

For research use only .