

产品名称: Cyclophilin B Mouse Monoclonal Antibody
产品货号: AMM85037



产品概述 (Summary)

产品名称 (Production Name)	Cyclophilin B Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.5%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Cyclophilin B
别名 (Alternative Names)	PPIB; CYPB; Peptidyl-prolyl cis-trans isomerase B; PPIase B; CYP-S1; Cyclophilin B; Rotamase B; S-cyclophilin; SCYLP
基因 ID (Gene ID)	5479.0
蛋白 ID (SwissProt ID)	P23284.Synthetic peptide conjugated to KLH.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 24 kDa; Observed MW: 24 kDa

研究背景 (Background)


产品名称: Cyclophilin B Mouse Monoclonal Antibody
产品货号: AMM85037




PPlases accelerate the folding of proteins. It catalyzes the cis-trans isomerization of proline imidic peptide bonds in oligopeptides.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of Cyclophilin B in PC-3, T47D, COS7, C6 lysates using Cyclophilin B antibody




Western blot analysis of Cyclophilin B in HeLa, COS7, PC-3, Human placenta, Jurkat lysates using Cyclophilin B antibody.



Immunohistochemistry analysis of paraffin-embedded Human Hepatocarcinoma using Cyclophilin B antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: Cyclophilin B Mouse Monoclonal Antibody
产品货号: AMM85037



Immunohistochemical analysis of paraffin-embedded Human tonsils using Cyclophilin B antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .