

产品名称: CHD3 Mouse Monoclonal Antibody
产品货号: AMM85042



产品概述 (Summary)

产品名称 (Production Name)	CHD3 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB, ICC, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CHD3 CHD3; Chromodomain-helicase-DNA-binding protein 3; CHD-3; ATP-dependent helicase CHD3; Mi-2 autoantigen 240 kDa protein; Mi2-alpha; Zinc finger helicase; hZFH
别名 (Alternative Names)	
基因 ID (Gene ID)	1107.0
蛋白 ID (SwissProt ID)	Q12873. Purified recombinant human CHD3 (C-terminus) protein fragments expressed in E.coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:50-1:200, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 227 kDa; Observed MW: 260 kDa

产品名称: CHD3 Mouse Monoclonal Antibody
产品货号: AMM85042




研究背景 (Background)


Component of the histone deacetylase NuRD complex which participates in the remodeling of chromatin by deacetylating histones. Required for anchoring centrosomal pericentrin in both interphase and mitosis, for spindle organization and centrosome integrity.

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of CHD3 (Cterminus) in NTERA-2, 3T3, C6, K562 and HeLa lysates using CHD3 (Cterminus) antibody.



Immunoprecipitation analysis of CHD3 in HeLa lysates using CHD3 antibody.

产品名称: CHD3 Mouse Monoclonal Antibody
产品货号: AMM85042



Immunocytochemistry analysis of CHD3 in HeLa using CHD3 (Cterminus) antibody.

注意事项 (Note)

For research use only .