

产品名称: ICAM1 Mouse Monoclonal Antibody
产品货号: AMM85044



产品概述 (Summary)

产品名称 (Production Name)	ICAM1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ICAM1
别名 (Alternative Names)	ICAM1; Intercellular adhesion molecule 1; ICAM-1; Major group rhinovirus receptor; CD antigen CD54
基因 ID (Gene ID)	3383.0
蛋白 ID (SwissProt ID)	P05362. Purified recombinant human CD54(ICAM-1) protein fragments expressed in E.coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 58 kDa; Observed MW: 96 kDa

产品名称: ICAM1 Mouse Monoclonal Antibody
产品货号: AMM85044




研究背景 (Background)


ICAM1 a type I membrane protein of the immunoglobulin superfamily. Is a ligand for the leukocyte adhesion LFA-1 protein (Integrin alpha-L/beta-2) and a Rhinovirus receptor. Typically expressed on endothelial cells and cells of the immune system.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of CD54(ICAM1) in SW480 and Raji lysates using CD54(ICAM1) antibody.



Immunocytochemistry analysis of ICAM1 in SW480 using CD54(ICAM1) antibody.

注意事项 (Note)

For research use only .