

产品名称: FOXP3 Mouse Monoclonal Antibody
产品货号: AMM85047



产品概述 (Summary)

产品名称 (Production Name)	FOXP3 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Transfected

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FOXP3
别名 (Alternative Names)	FOXP3; IPEX; JM2; Forkhead box protein P3; Scurfin
基因 ID (Gene ID)	50943.0
蛋白 ID (SwissProt ID)	Q9BZS1. Purified recombinant human FOXP3 protein fragments expressed in E.coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 47 kDa; Observed MW: 47 kDa

研究背景 (Background)


产品名称: FOXP3 Mouse Monoclonal Antibody
产品货号: AMM85047




Defects in FOXP3 are the cause of immunodeficiency polyendocrinopathy, enteropathy, X-linked syndrome (IPEX) [MIM:304790]; also known as X-linked autoimmunity-immunodeficiency syndrome. IPEX is characterized by neonatal onset insulin-dependent diabetes mellitus, infections, secretory diarrhea, thrombocytopenia, anemia and eczema. It is usually lethal in infancy.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of FOXP3 in CHO-K1 lysates(A) and CHO-K1 transfected by FOXP3(B) lysates using FOXP3 antibody.



Immunocytochemistry analysis of FOXP3 mouse mAbCHO-K1 which transfected by FOXP3 using FOXP3 antibody.

注意事项 (Note)

For research use only .