

产品名称: GRB2 Mouse Monoclonal Antibody
产品货号: AMM85050



产品概述 (Summary)

产品名称 (Production Name)	GRB2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.5%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GRB2
别名 (Alternative Names)	GRB2; ASH; Growth factor receptor-bound protein 2; Adapter protein GRB2; Protein Ash; SH2/SH3 adapter GRB2
基因 ID (Gene ID)	2885.0
蛋白 ID (SwissProt ID)	P62993.Recombinant protein of human GRB2.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 25 kDa; Observed MW: 25 kDa

研究背景 (Background)

产品名称: GRB2 Mouse Monoclonal Antibody
产品货号: AMM85050




Adapter protein that provides a critical link between cell surface growth factor receptors and the Ras signaling pathway.

研究领域 (Research Area)

TGF-beta signaling pathway, PI3K-Akt signaling pathway, MAPK signaling pathway

图片 (Image Data)



Western blot analysis of GRB2 in K562, Raji, Jurkat, Ramos, C6, 3T3 and HeLa lysates using GRB2 antibody.

注意事项 (Note)

For research use only .