

产品名称: CDX2 Mouse Monoclonal Antibody
产品货号: AMM85057



产品概述 (Summary)

产品名称 (Production Name)	CDX2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB, ICC, IP
种属反应性 (Reactivity)	Human, Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CDX2
别名 (Alternative Names)	CDX2; CDX3; Homeobox protein CDX-2; CDX-3; Caudal-type homeobox protein 2
基因 ID (Gene ID)	1045.0
蛋白 ID (SwissProt ID)	Q99626. Purified recombinant human CDX2 protein fragments expressed in E.coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:50-1:200, IP 1:10-1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 34 kDa; Observed MW: 34 kDa

产品名称: CDX2 Mouse Monoclonal Antibody
产品货号: AMM85057




研究背景 (Background)


Involved in the transcriptional regulation of multiple genes expressed in the intestinal epithelium. Important in broad range of functions from early differentiation to maintenance of the intestinal epithelial lining of both the small and large intestine.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of CDX2 in HeLa, COS7 and K562 lysates using CDX2 antibody.



Immunoprecipitation analysis of CDX2 in HeLa lysates using CDX2 antibody.



Immunocytochemistry analysis of CDX2 in HeLa using CDX2 antibody.

产品名称: CDX2 Mouse Monoclonal Antibody
产品货号: AMM85057



注意事项 (Note)

For research use only .