

产品名称: 14-3-3 eta Mouse Monoclonal Antibody
产品货号: AMM85062



产品概述 (Summary)

产品名称 (Production Name)	14-3-3 eta Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	14-3-3 eta 14-3-3 protein eta; 1433F_HUMAN; Brain protein 14-3-3; eta isoform;
别名 (Alternative Names)	HGNC:12853; Protein AS1; Tyrosine 3/tryptophan 5 monooxygenase activation protein eta polypeptide; YWHA 1; YWHA1; Ywhah.
基因 ID (Gene ID)	7533.0
蛋白 ID (SwissProt ID)	Q04917. Antibodies are produced by immunizing animals with a synthetic peptide corresponding to the sequence of human 14-3-3 η.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 28 kDa; Observed MW: 28 kDa

产品名称: 14-3-3 eta Mouse Monoclonal Antibody
产品货号: AMM85062




研究背景 (Background)

Adapter protein implicated in the regulation of a large spectrum of both general and specialized signaling pathways. Binds to a large number of partners, usually by recognition of a phosphoserine or phosphothreonine motif. Binding generally results in the modulation of the activity of the binding partner.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of 14-3-3 eta/YWHAH in HeLa, 3T3, COS7 and C2C12 lysates using 14-3-3 eta/YWHAH antibody.

注意事项 (Note)

For research use only .