

产品名称: Lyn Mouse Monoclonal Antibody
产品货号: AMM85070



产品概述 (Summary)

产品名称 (Production Name)	Lyn Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Lyn LYN; JTK8; Tyrosine-protein kinase Lyn; Lck/Yes-related novel protein tyrosine kinase; V-yes-1 Yamaguchi sarcoma viral related oncogene homolog; p53Lyn; p56Lyn
别名 (Alternative Names)	
基因 ID (Gene ID)	4067.0
蛋白 ID (SwissProt ID)	P07948. Purified recombinant human Lyn protein fragments expressed in E.coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 59 kDa; Observed MW: 56,53 kDa

产品名称: Lyn Mouse Monoclonal Antibody
产品货号: AMM85070




研究背景 (Background)

Lyn a tyrosine kinase of the Src family that plays a role in B cell antigen-receptor signaling.

研究领域 (Research Area)

Jak-STAT signaling pathway

图片 (Image Data)



Western blot analysis of Lyn in Ramos, Raji, K562 and HeLa lysates using Lyn antibody.

注意事项 (Note)

For research use only .