

产品名称: NSE Mouse Monoclonal Antibody
产品货号: AMM85072



产品概述 (Summary)

产品名称 (Production Name)	NSE Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ICC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.5%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NSE
别名 (Alternative Names)	ENO2; Gamma-enolase; 2-phospho-D-glycerate hydro-lyase; Enolase 2; Neural enolase; Neuron-specific enolase; NSE
基因 ID (Gene ID)	2026.0
蛋白 ID (SwissProt ID)	P09104.Synthetic Peptide of NSE

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 47 kDa; Observed MW: 47 kDa

研究背景 (Background)


产品名称: NSE Mouse Monoclonal Antibody
产品货号: AMM85072




ENO2 an enzyme with 2-phospho-D-glycerate hydro-lyase activity. One of the three enolase isoenzymes found in mammals. This isoenzyme, a homodimer, is found in mature neurons and cells of neuronal origin. A switch from alpha enolase to gamma enolase occurs in neural tissue during development in rats and primates.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of Enolase2 in COS7, C6 and HeLa lysates using Enolase2 antibody.




Western blot analysis of NSE in Enolase in CHO-K1, rat brain, C6, 3T3, HeLa lysates using Enolase antibody.



Immunocytochemistry analysis of NSE (green) in HeLa using NSE antibody ,and DAPI(blue)

产品名称: NSE Mouse Monoclonal Antibody
产品货号: AMM85072



Immunohistochemistry analysis of paraffin-embedded Human small cell carcinoma of lung tissue using NSE antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .