

产品名称: FH Mouse Monoclonal Antibody
产品货号: AMM85073



产品概述 (Summary)

产品名称 (Production Name)	FH Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ICC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.5%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FH
别名 (Alternative Names)	Fumarate hydratase; mitochondrial; Fumarase
基因 ID (Gene ID)	2271.0
蛋白 ID (SwissProt ID)	P07954.Synthetic Peptide of FH

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 55 kDa; Observed MW: 50 kDa

研究背景 (Background)

Also acts as a tumor suppressor.MiscellaneousThere are 2 substrate-binding sites: the catalytic A site, and the non-catalytic


产品名称: FH Mouse Monoclonal Antibody
产品货号: AMM85073




B site that may play a role in the transfer of substrate or product between the active site and the solvent. Alternatively, the B site may bind allosteric effectors .

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of Fumarase in MCF-7, Jurkat and HeLa lysates using Fumarase antibody.




Western blot analysis of FH in K562, Ramos, Jurkat, HeLa lysates using Fumarase antibody.



Immunocytochemistry analysis of FH/Fumarase (green) in HeLa using FH/Fumarase antibody ,and DAPI(blue)

产品名称: FH Mouse Monoclonal Antibody
产品货号: AMM85073



Immunohistochemistry analysis of paraffin-embedded Human tonsil tissue using FH/Fumarase antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .