

产品名称: ABCB5 Mouse Monoclonal Antibody
产品货号: AMM85074



产品概述 (Summary)

产品名称 (Production Name)	ABCB5 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ABCB5
别名 (Alternative Names)	ABCB5 P-gp; ATP binding cassette sub family B (MDR/TAP) member 5; P glycoprotein ABCB5
基因 ID (Gene ID)	340273.0
蛋白 ID (SwissProt ID)	Q2M3G0.Synthetic Peptide of ABCB5

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 139 kDa; Observed MW: 90 kDa

研究背景 (Background)


产品名称: ABCB5 Mouse Monoclonal Antibody
产品货号: AMM85074




Drug efflux transporter present in a number of stem cells that acts as a regulator of cellular differentiation. Able to mediate efflux from cells of the rhodamine dye and of the therapeutic drug doxorubicin. Specifically present in limbal stem cells, where it plays a key role in corneal development and repair.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of ABCB5 in K562 lysates using ABCB5 antibody.



Western blot analysis of ABCB5 in HepG2, K562 lysates using ABCB5 antibody.



Immunocytochemistry analysis of ABCB5 (green) in Hela using ABCB5 antibody ,and DAPI(blue)

产品名称: ABCB5 Mouse Monoclonal Antibody
产品货号: AMM85074



注意事项 (Note)

For research use only .