

产品名称: Bcl2 Mouse Monoclonal Antibody
产品货号: AMM85085



产品概述 (Summary)

产品名称 (Production Name)	Bcl2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ICC
种属反应性 (Reactivity)	Human,Other

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.5%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Bcl2
别名 (Alternative Names)	BCL2; Apoptosis regulator Bcl-2
基因 ID (Gene ID)	596.0
蛋白 ID (SwissProt ID)	P10415.Synthetic Peptide of Bcl-2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 26 kDa; Observed MW: 26 kDa

研究背景 (Background)

This gene encodes an integral outer mitochondrial membrane protein that blocks the apoptotic death of some cells such as

产品名称: Bcl2 Mouse Monoclonal Antibody
产品货号: AMM85085




lymphocytes. Constitutive expression of BCL2, such as in the case of translocation of BCL2 to Ig heavy chain locus, is thought to be the cause of follicular lymphoma.


研究领域 (Research Area)

Apoptosis, TGF-beta signaling pathway, PI3K-Akt signaling pathway

图片 (Image Data)



Western blot analysis of Bcl2 in chicken lysates using Bcl2 antibody




Western blot analysis of Bcl2 in Human breast cancer cell line MCF-7(A), MDAMB231(B) and Cal51(C) using Bcl2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil tissue using Bcl2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: Bcl2 Mouse Monoclonal Antibody
产品货号: AMM85085



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Bcl2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .