

产品名称: PRMT6 Mouse Monoclonal Antibody
产品货号: AMM85087



产品概述 (Summary)

产品名称 (Production Name)	PRMT6 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ICC
种属反应性 (Reactivity)	Human,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.5%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PRMT6
别名 (Alternative Names)	HRMT1L6; PRMT 6
基因 ID (Gene ID)	55170.0
蛋白 ID (SwissProt ID)	Q96LA8.Purified recombinant human PRMT6 protein fragments expressed in E.coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 42 kDa; Observed MW: 42 kDa

研究背景 (Background)


产品名称: PRMT6 Mouse Monoclonal Antibody
产品货号: AMM85087



Arginine methyltransferase that can catalyze the formation of both omega-N monomethylarginine (MMA) and asymmetrical dimethylarginine (aDMA), with a strong preference for the formation of aDMA. Preferentially methylates arginyl residues present in a glycine and arginine-rich domain and displays preference for monomethylated substrates.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of PRMT6 in COS7, Jurkat, C6, 3T3 and HeLa lysates using PRMT6 antibody.

注意事项 (Note)

For research use only .