

产品名称: ATIC Mouse Monoclonal Antibody
产品货号: AMM85091



产品概述 (Summary)

产品名称 (Production Name)	ATIC Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Monkey,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.5%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ATIC
别名 (Alternative Names)	PURH; AICAR; AICARFT; IMPCHASE
基因 ID (Gene ID)	471.0
蛋白 ID (SwissProt ID)	P31939.Purified recombinant human ATIC protein fragments expressed in E.coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 65 kDa; Observed MW: 65 kDa

研究背景 (Background)


产品名称: ATIC Mouse Monoclonal Antibody
产品货号: AMM85091



Bifunctional enzyme that catalyzes 2 steps in purine biosynthesis.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of ATIC in Jurkat, C6, 3T3, COS7 and HeLa lysates using ATIC antibody.

注意事项 (Note)

For research use only .