

产品名称: DDX3 Mouse Monoclonal Antibody
产品货号: AMM85097



产品概述 (Summary)

产品名称 (Production Name)	DDX3 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB, ICC, IP
种属反应性 (Reactivity)	Human, Mouse, Rat, Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DDX3
别名 (Alternative Names)	DBX; DDX3; HLP2; DDX14.
基因 ID (Gene ID)	1654.0
蛋白 ID (SwissProt ID)	O00571. Purified recombinant human DDX3 protein fragments expressed in E.coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 73 kDa; Observed MW: 73 kDa

研究背景 (Background)

产品名称: DDX3 Mouse Monoclonal Antibody
产品货号: AMM85097




Multifunctional ATP-dependent RNA helicase.

研究领域 (Research Area)

Apoptosis

图片 (Image Data)



Western blot analysis of DDX3 in COS7, Jurkat, C6, 3T3 and HeLa lysates using DDX3 antibody.

注意事项 (Note)

For research use only .