

产品名称: H1L Mouse Monoclonal Antibody
产品货号: AMM85906



产品概述 (Summary)

产品名称 (Production Name)	H1L Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	VV、Vaccinia virus

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	H1L
别名 (Alternative Names)	Dual specificity protein phosphatase H1, 313-, Late protein H1, H1L
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P20495. Purified recombinant H1L fusion protein was used to produced this monoclonal antibody.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	19.7kDa

研究背景 (Background)


This Protein shows activity toward tyrosine-protein phosphatase as well as with serine-protein phosphatase.

产品名称: H1L Mouse Monoclonal Antibody
产品货号: AMM85906



研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of anti-H1L Mab in recombinant pGEX-H1L protein. H1L(arrow) was detected using the purified Mab.

注意事项 (Note)

For research use only .