

产品名称: CD34 Mouse Monoclonal Antibody
产品货号: AMM85908



产品概述 (Summary)

产品名称 (Production Name)	CD34 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD34
别名 (Alternative Names)	Hematopoietic progenitor cell antigen CD34, CD34, CD34
基因 ID (Gene ID)	947.0
蛋白 ID (SwissProt ID)	P28906. This CD34 antibody is generated from mice immunized with a KLH conjugated synthetic peptide between 185-214 amino acids from the Central region of human CD34.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, IHC 1:500-1:2000
蛋白分子量 (Molecular Weight)	40.7kDa

研究背景 (Background)


产品名称: CD34 Mouse Monoclonal Antibody
产品货号: AMM85908



CD34 is a monomeric cell surface antigen with a molecular mass of approximately 110 kD that is selectively expressed on human hematopoietic progenitor cells.

研究领域 (Research Area)

图片 (Image Data)



Anti-CD34 Antibody (Center) at 1:1000 dilution + NIH/3T3 whole cell lysate

注意事项 (Note)

For research use only .