

产品名称: MAPK15 Mouse Monoclonal Antibody
产品货号: AMM85913



产品概述 (Summary)

产品名称 (Production Name)	MAPK15 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MAPK15
别名 (Alternative Names)	Mitogen-activated protein kinase 15, MAP kinase 15, MAPK 15, Extracellular signal-regulated kinase 7, ERK-7, Extracellular signal-regulated kinase 8, ERK-8, MAPK15, ERK7, ERK8
基因 ID (Gene ID)	225689.0
蛋白 ID (SwissProt ID)	Q8TD08. This MAPK15 monoclonal antibody is generated from mouse immunized with MAPK15 recombinant protein.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC 1:20-1:50
蛋白分子量 (Molecular Weight)	59.8kDa

产品名称: MAPK15 Mouse Monoclonal Antibody
产品货号: AMM85913




研究背景 (Background)

In vitro, phosphorylates MBP.


研究领域 (Research Area)

Jak-STAT signaling pathway

图片 (Image Data)



MAPK15 Antibody (Cat. #AMM85913) western blot analysis in MCF-7 cell line lysates (35µg/lane). This demonstrates the MAPK15 antibody detected the MAPK15 protein (arrow).



Confocal immunofluorescent analysis of MAPK15 Antibody (Cat#AMM85913) with MCF-7 cell followed by Alexa Fluor® 488-conjugated goat anti-mouse IgG (green). DAPI was used to stain the cell nuclear (blue).

注意事项 (Note)

For research use only .