

产品名称: MIPEP Mouse Monoclonal Antibody
产品货号: AMM85914



产品概述 (Summary)

产品名称 (Production Name)	MIPEP Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MIPEP
别名 (Alternative Names)	Mitochondrial intermediate peptidase, MIP, MIPEP, MIP
基因 ID (Gene ID)	4285.0
蛋白 ID (SwissProt ID)	Q99797. This MIPEP monoclonal antibody is generated from mouse immunized with MIPEP recombinant protein.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:200-1:1000
蛋白分子量 (Molecular Weight)	80.4kDa

研究背景 (Background)


Cleaves proteins, imported into the mitochondrion, to their mature size.

产品名称: MIPEP Mouse Monoclonal Antibody
产品货号: AMM85914



研究领域 (Research Area)

图片 (Image Data)



MIPEP Antibody (Cat. #AMM85914) western blot analysis in 293 cell line lysates (35 μ g/lane). This demonstrates the MIPEP antibody detected the MIPEP protein (arrow).

注意事项 (Note)

For research use only .