

产品名称: MASTL Mouse Monoclonal Antibody
产品货号: AMM85918



产品概述 (Summary)

产品名称 (Production Name)	MASTL Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MASTL
别名 (Alternative Names)	Serine/threonine-protein kinase greatwall, GW, GWL, hGWL, Microtubule-associated serine/threonine-protein kinase-like, MAST-L, MASTL, GW, GWL, THC2
基因 ID (Gene ID)	84930.0
蛋白 ID (SwissProt ID)	Q96GX5. This MASTL monoclonal antibody is generated from mouse immunized with MASTL recombinant protein.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:15000
蛋白分子量 (Molecular Weight)	97.3kDa

产品名称: MASTL Mouse Monoclonal Antibody
产品货号: AMM85918




研究背景 (Background)

This gene encodes a microtubule-associated serine/threonine kinase. Mutations at this locus have been associated with autosomal dominant thrombocytopenia, also known as thrombocytopenia-2. Alternatively spliced transcript variants have been described for this locus.

研究领域 (Research Area)

图片 (Image Data)



MASTL (Cat. #AMM85918) western blot analysis in HepG2 cell line lysates (35µg/lane). This demonstrates the MASTL antibody detected the MASTL protein (arrow).

注意事项 (Note)

For research use only .