

产品名称: Rat Serpinh1 Mouse Monoclonal Antibody
产品货号: AMM85919



产品概述 (Summary)

产品名称 (Production Name)	Rat Serpinh1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Rat Serpinh1
别名 (Alternative Names)	Serpin H1, 47 kDa heat shock protein, Collagen-binding protein, Colligin, GP46, Serpinh1, Cbp1, Hsp47
基因 ID (Gene ID)	29345.0
蛋白 ID (SwissProt ID)	P29457. This Rat Serpinh1 antibody is generated from mice immunized with a KLH conjugated synthetic peptide between 222-251 amino acids from the Central region of rat Serpinh1.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:4000, IHC 1:100-1:500
蛋白分子量 (Molecular Weight)	46.5kDa


产品名称: Rat Serpinh1 Mouse Monoclonal Antibody
产品货号: AMM85919



研究背景 (Background)

研究领域 (Research Area)

图片 (Image Data)



Anti-Rat Serpinh1 Antibody (Center) (Ascites) at 1:4000 dilution + rat heart lysate

注意事项 (Note)

For research use only .