

产品名称: LHX2 Mouse Monoclonal Antibody  
产品货号: AMM85920



## 产品概述 (Summary)

|                        |                                |
|------------------------|--------------------------------|
| 产品名称 (Production Name) | LHX2 Mouse Monoclonal Antibody |
| 描述 (Description)       | Mouse monoclonal Antibody      |
| 宿主 (Host)              | Mouse                          |
| 应用 (Application)       | WB                             |
| 种属反应性 (Reactivity)     | Human                          |

## 产品性能 (Performance)

|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                             |
| 修饰 (Modification)   | Unmodified                                                                               |
| 同种型 (Isotype)       | Mouse IgM                                                                                |
| 克隆 (Clonality)      | Monoclonal                                                                               |
| 形式 (Form)           | Liquid                                                                                   |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| 储存溶液 (Buffer)       | Purified antibody in PBS with 0.05% sodium azide.                                        |
| 纯化方式 (Purification) | Affinity Purification                                                                    |

## 免疫原信息 (Immunogen)

|                        |                                                                                                        |
|------------------------|--------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | LHX2                                                                                                   |
| 别名 (Alternative Names) | LIM/homeobox protein Lhx2, Homeobox protein LH-2, LIM homeobox protein 2, LHX2, LH2                    |
| 基因 ID (Gene ID)        | 9355.0                                                                                                 |
| 蛋白 ID (SwissProt ID)   | P50458. This LHX2 monoclonal antibody is generated from mouse immunized with LHX2 recombinant protein. |

## 产品应用 (Application)

|                          |                 |
|--------------------------|-----------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:1000 |
| 蛋白分子量 (Molecular Weight) | 44.4kDa         |

## 研究背景 (Background)


**产品名称: LHX2 Mouse Monoclonal Antibody**  
**产品货号: AMM85920**



This gene encodes a protein belonging to a large protein family, members of which carry the LIM domain, a unique cysteine-rich zinc-binding domain. The encoded protein may function as a transcriptional regulator. The protein can recapitulate or rescue phenotypes in *Drosophila* caused by a related protein, suggesting conservation of function during evolution. [provided by RefSeq].

## 研究领域 (Research Area)

## 图片 (Image Data)



LHX2 Antibody (Cat. #AMM85920) western blot analysis in MDA-MB435 cell line lysates (35µg/lane). This demonstrates the LHX2 antibody detected the LHX2 protein (arrow).

## 注意事项 (Note)

For research use only .