

产品名称: GSN (N-term) Mouse Monoclonal Antibody
产品货号: AMM85921



产品概述 (Summary)

产品名称 (Production Name)	GSN (N-term) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgM
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GSN (N-term)
别名 (Alternative Names)	Gelsolin, AGEL, Actin-depolymerizing factor, ADF, Brevin, GSN
基因 ID (Gene ID)	2934.0
蛋白 ID (SwissProt ID)	P06396. This GSN antibody is generated from mice immunized with a KLH conjugated synthetic peptide between 230-259 amino acids from the N-terminal region of human GSN.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:200-1:4000
蛋白分子量 (Molecular Weight)	85.7kDa

研究背景 (Background)


产品名称: GSN (N-term) Mouse Monoclonal Antibody
产品货号: AMM85921



The protein encoded by this gene binds to the 'plus' ends of actin monomers and filaments to prevent monomer exchange. The encoded calcium-regulated protein functions in both assembly and disassembly of actin filaments. Defects in this gene are a cause of familial amyloidosis Finnish type (FAF). Multiple transcript variants encoding several different isoforms have been found for this gene.

研究领域 (Research Area)

图片 (Image Data)



GSN Antibody (N-term) (Cat. #AMM85921) western blot analysis in MDA-MB435 cell line lysates (35µg/lane). This demonstrates the GSN antibody detected the GSN protein (arrow).

注意事项 (Note)

For research use only .